




Agilent OpenLAB Chromatography Data System (CDS)

EZChrom and ChemStation Editions



Hardware and Software Requirements



Agilent Technologies

Notices

© Agilent Technologies, Inc. 2011

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

M8301-90041

Edition

Second edition, March 2011

Printed in USA

Agilent Technologies, Inc.

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as "Commercial computer software" as defined in DFAR 252.227-7014 (June 1995), or as a "commercial item" as defined in FAR 2.101(a) or as "Restricted computer software" as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies' standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Contents

Introduction 1

Hardware 1

Memory, processing, disk and communications 1

Agilent Instrument Controller (AIC) compatibility 2

Number of Instruments Supported 2

Software 3

Operating systems 3

OpenLAB ECM software 3

Licensing 4

Application virtualization/delivery 4

Other software 4

Language compatibility 5

Important notes 5

Introduction

This document details the minimum hardware and software requirements for supporting the Agilent OpenLAB Chromatography Data System (CDS) family of products. This document is valid for OpenLAB CDS Shared Services A.01.02, ChemStation Edition C.01.02, and EZChrom Edition A.04.02.

Hardware

Memory, processing, disk and communications

Table 1 lists minimum and recommended hardware configurations by category and OpenLAB CDS software product.

Table 1 Minimum and recommended hardware configurations

Item	Workstation, Client, or EZChrom AIC on 32bit OS	Workstation, Client, or EZChrom AIC on Win7 or 64bit OS	ChemStation AIC on 64bit OS	Shared Services Server (32 or 64 bit OS)
Processor Speed (CPU)				
Minimum	2 GHz	2 GHz	2GHz	2 GHz
Recommended	3 GHz Dual core	3 GHz Dual core	Dual Quad Core	3 GHz Dual core
Physical Memory (RAM)				
Minimum	1 GB	2 GB	8GB	4 GB
Recommended	4 GB	4GB (32-bit), 6GB (64-bit)	10 GB	6GB (64-bit)
Disk Space Minimum				
Total	80 GB	80 GB	250 GB	80 GB
Available	40GB	40GB	50 GB	40GB
Optical Drive	DVD ROM (12x)	DVD ROM (12x)	DVD ROM (12x)	DVD ROM (12x)
Network Adapter				
Minimum	10/100 LAN Adapter	10/100 LAN Adapter	10/100 LAN Adapter	10/100 LAN Adapter
Recommended	100/1000 LAN Adapter	100/1000 LAN Adapter	100/1000 LAN Adapter	100/1000 LAN Adapter

Agilent Instrument Controller (AIC) compatibility

Table 2 and **Table 3** list the Agilent Instrument Controllers that support installation of OpenLAB CDS EZChrom Edition or ChemStation Edition.

Table 2 Agilent Instrument Controllers that support installation of OpenLAB CDS EZChrom Edition

AIC Model	CPU	RAM	Disk Space
5000	2.8 GHz	1 GB	80 GB
7800	2.3 GHz Dual-Core	2 GB	135 GB
7900	3.0 GHz Dual-Core	2 GB	232 GB (2x116)
8000	3.16 GHz Dual-Core	4 GB	250 GB

Revision A and B Agilent instrument controllers are not supported for use with Agilent OpenLAB CDS.

Table 3 Agilent Instrument Controllers that support installation of OpenLAB CDS ChemStation Edition

AIC Model	CPU	RAM	Disk Space
ML150 (P/N M8304AA)	2.00 GHz 2xQuad Core	8 GB	2x500 GB

Number of Instruments Supported

There is a limit to the number of instruments (for example, GCs, LCs or CEs) that can be configured on a single OpenLAB CDS **ChemStation Edition** WorkStation or Agilent Instrument Controller (AIC), as shown in **Table 4**.

Each instrument can consist of several modules, as in the case of Agilent modular LC systems. The theoretical maximum number of modules is 31; the recommended maximum number is 14.

Table 4 Maximum Number of Instruments: ChemStation Edition

	WorkStation	AIC
GC, LC (2D)	4	10
LC (3D), CE	2	5
LCMS, CEMS	1	Not Supported

There is no limit to the number of AICs in a distributed system.

Similarly, **Table 5** below shows the instrument connection limit for workstations and AICs running OpenLAB CDS **EZChrom Edition**.

Table 5 Maximum Number of Instruments: EZChrom Edition

	WorkStation	AIC
GC, LC (2D)	4	4
LC (3D), CE	2	2
LCMS, CEMS	Not Supported	Not Supported

Software

The following information applies to **all** installable OpenLAB CDS components except where otherwise noted.

Operating systems

The following Operating systems are supported for the installation and use of the specified OpenLAB CDS components:

Component	Windows XP SP3	Windows Vista SP1	Windows 7 (32/64bit)	Windows Server 2003 R2 SP2	Windows Server 2008 R2
Shared Services Server	✗	✗	✗	✓	✓
ChemStation Edition Workstation/NWS	✓	✓	✓	✗	✗
ChemStation Edition AIC	✗	✗	✗	✗	✓
EZChrom Edition Workstation/NWS	✓	✓	✓	✗	✗
EZChrom Edition Client	✓	✓	✓	✓	✓
EZChrom Edition AIC	✓	✓	✓	✗	✗

OpenLAB ECM software

In order for OpenLAB CDS to use ECM as a file storage location, you will need these ECM product revisions:

- Agilent OpenLAB ECM Server 3.3.2 SP1-3.4.1
- Agilent OpenLAB ECM API v. 1.0.25.x (Available as ECM 3.3.2 Hotfix03, included with ECM 3.4.1)

Licensing

Agilent OpenLAB CDS supports the use of a central licensing server for the distribution and tracking of license entitlements. The following software is supported for this purpose:

- Flex-Net Publisher v. 11.8

This software is installed with the installation of OpenLAB Shared Services Server components.

Databases

Agilent OpenLAB CDS Shared Services Servers manage information using a database. The database is installed and configured automatically during OpenLAB CDS workstation installation.

When deploying networked system configurations, the database requires additional installation to ensure acceptable performance of the system.

The following database software is supported for hosting shared services databases:

- SQL Server 2008 R2 Express (32-bit) : Installed with product
- SQL Server 2008 R2 Standard or Enterprise (32-bit)
- SQL Server 2008 R2 Standard or Enterprise (64-bit) (Windows 2008 Server R2 databases only)
- Oracle 11g Standard and Enterprise Editions

Note: It is strongly recommended that you store database files and transaction logs (when applicable) on physically exclusive hard drives.

Application virtualization/delivery

OpenLAB CDS shared services servers, instrument controllers and client systems can be run on virtual machines that meet the requirements outlined above. Use the VMWare vSphere 4.0 software for host machines. See Network Requirements guides for restrictions regarding implementation of virtualized instrument controllers.

OpenLAB CDS EZChrom Edition clients are also supported to run via terminal services and Citrix XenApp 5.0.

Other software

The following software must be installed on any supported operating system prior to installing OpenLAB CDS components:

- Microsoft .Net Framework 3.5 with Service Pack 1
- Windows Installer 4.5 (required for SQL Installations on Shared Servers)

Language compatibility

The English version of OpenLAB CDS is validated on Windows English and Western European language versions.

Localized versions of OpenLAB CDS are supported on localized language versions of Windows, using default system fonts (Chinese: Simsun; Japanese: MS UI Gothic; Portuguese: MS Sans Serif).

Non-localized instrument drivers are supported, and will appear in English even when running localized versions of OpenLAB CDS.

Important notes

- Unless otherwise stated, workstations, clients and Agilent instrument controllers listed above are tested and supported with any combination of 4 instruments (LC's containing PDA/DAD count as 2 instruments) with up to one Headspace or CTC Sampler running simultaneously.
- Disk space requirements should be adjusted based on the **number and type of instruments** and **archival periodicity**. Agilent recommends a period of expected disk usage one year beyond the recommendation of the operating system.

Typical file sizes for reference:

- 2D data: 60min, 10Hz, 2 channel data ~600Kb
- 3D data: 60min, 10Hz, 5 channel plus spectra at 1nm resolution from 200-400nm ~300Mb
- Citrix and Terminal Servers are supported for running OpenLAB CDS EZChrom Edition clients only.